

GENERAL INFORMATION FOR SAMPLE SUBMISSION

Necessary documents and regulatory essentials



SHIPPING FROM EU MEMBER STATES

Submitting original/native samples (e. g. swabs, serum, tissue, feces etc.)

Samples must be accompanied by a commercial document and by our submission form.

The CN-code on the commercial document can be indicated as follows:

CN-code for isolates: 30024900

CN-code for blood samples: 30029030

CN-code for other diagnostic material: 05119985

Submitting FTA-cards/AniCards

Cards should be accompanied by our submission form.

 **Submission form incl. commercial document and packaging instructions**



PACKAGING

Samples sent to our laboratory should only contain Category B material acc. to UN3373 which includes – among others – animal material that gets transported for the purpose of diagnosis or examination. Such material includes excreta and blood as well as fluids and other tissues.

Original/native samples belonging to the Category B should be packed according to the packaging instruction P650, which implies to UN 3373. The package should consist of three components:

- a primary leakage proof packaging
- a secondary packaging with sufficient absorbing material
- an outer packing

For further information regarding the packaging see the packaging instructions on our website:



www.anicon.eu/swine-downloads



www.anicon.eu/poultry-downloads



SHIPPING FROM NON-EU MEMBER STATES

Submitting original/native samples (e. g. swabs, serum, tissue, feces etc.)

Sending original samples from non-EU-countries is restricted. Please contact us in any case for possibilities and further information.

Submitting FTA-cards/AniCards

Cards should be accompanied by our submission form.

For the shipment of FTA-cards/AniCards from non-EU member states, please indicate the following HS-code and EORI-number to the shipping-company:

HS-code: 48119000

EORI-number: DE6173438

For further information and support please do not hesitate to contact us. We are at your disposal!

E-mail: office-de@san-group.com · Phone: +49 4473 9438 30

Muehlenstraße 13 | 49685 Hoeltinghausen

ANICON®