

Veterinary Sample Submission Card



INSTRUCTIONS:

- 1 Label the AniCard with sample ID.
- 2 Note sample description.
- 3 Apply one sample to each spot.
- 4 Let AniCard air-dry.
- 5 Place AniCard into the bag and seal.
- 6 Fill in sample submission form.



Send samples with completed form to:
AniCon Labor GmbH
Muehlenstr. 13
49685 Hoeltinghausen, Germany

**For veterinary diagnostic use only.
Not available or to be used in the USA.**




For further information:
· Safety and import regulations
· User guide and manuals
· Sample submission form
Visit: anicon.eu/card


..... 5.
..... 4.
..... 3.
..... 2.
..... 1.
..... 6.
..... 7.
..... 8.
..... 9.
..... 10.

.....
Sample ID / Customer / Barcode


1




2




3




4



5




6




7



8



9



AniCon
AniCon Labor GmbH

10

