



Pathology

Selective sampling for bacteriology

Extended sampling efforts for PCR and / or cultural virology

Pathological-anatomical reporting

Coprosopic / microscopic examination of gut (mucosal scraping)

Parasitologic investigation (flotation/sedimentation)

Patho-histological examination (per tissue sample)

Material for sampling

FTA-cards (per card)

Reaction tube 1,5 ml (per 500 tubes)

Tube 13 ml and plug (per 500 tubes)

Salmonella-sampling swab (per 50 swabs)

Boot/Sock swab sampling material (per 100 swabs)

Sampling beaker sterile (per 5 beaker)

Swab for anaerobic bacteria (per 50 swabs)

Fine swab (aluminium) with transport medium (bacteriology) (per 100 swabs)

Fine swab (aluminium) without transport medium (bacteriology/PCR) (per 100 swabs)

Swab (wooden) without transport medium (bacteriology) (per 100 swabs)

Swab with preservative agent (Mykoplasma or Virus culture) (per 50 swabs)

Swab with preservative agent (Brachyspira) (per 50 swabs)

Paperstrips (ELISA) (per 100 strips)